

Gas Sampling Solutions

CO₂ CANNULA AND CO₂ MASK PROVIDE SOLUTIONS
FOR ROUTINE GAS MONITORING

Precision
Sampling



Westmed

When Your Care is Critical


Gas Sampling Solutions

CO₂ CANNULA AND CO₂ MASK
 PROVIDE SOLUTIONS FOR ROUTINE GAS MONITORING

The Westmed CO₂ Cannula, CO₂ Mask and CO₂ Accessories provide solutions to routine CO₂ monitoring while administering conscious sedation or during post operative recovery. Westmed offers a full array of products to meet your gas sampling requirements.

- Specific configurations to meet your requirements.
- Precision sampling of end-tidal CO₂ and anesthetic gases.
- Lightweight and economical
- Variety of lengths, ID's and luer fittings to accommodate a full range of monitors.
- Standard CO₂ lines for monitoring CO₂, O₂ and N₂O.
- Gas Sampling Line Accessories include Gas Sampling Elbow and Filter.

More Gas Sampling Products



Ask your local representative about other standard and custom configurations that are available to meet your gas monitoring needs.

Please log onto www.westmedinc.com for a complete listing of Westmed anesthesia and respiratory products.

© 2007 Westmed, Inc. All rights reserved. PN 74428-03

Part Number	Description	Case Quantity
Gas Sampling / Monitoring Lines		
0514	.040ID x 96" Male x Male	.50
0527	.040ID x 96" Male x Female	.50
0516	.047ID x 120" Male x Male Coextruded	.50
0517	.047ID x 120" Male x Female Coextruded	.50
0518	.040"OD X 120" Male x Male	.50
0519	.050ID x 72" Male x Male	.50
0521	.050ID x 72" Male x Female	.50
0512	.050ID x 120" Male x Male	.50
0524	.050ID x 120" Male x Male, 5.0 Micron Filter	.50
0522	.050ID x 120" Male x Male, .5 Micron Filter	.50
0513	.050ID x 120" Male x Female	.50
0525	.050ID x 120" Male x Female, 5.0 Micron Filter	.50
0523	.050ID x 120" Male x Female, .5 Micron Filter	.50
0510	.060ID x 120" Male x Female	.50
0511	.060ID x 120" Male x Male	.50
Gas Sampling Line Filters		
0529	0.5 Micron Sampling Line Filter, m/f Luer	.50
0520	0.8 Micron Sampling Line Filter m/f Luer	.100
0526	5.0 Micron Sampling Line Filter m/f Luer	.100
CO₂ Cannulas		
0503	CO ₂ Sampling Line (Male Luer)	25
0504	CO ₂ Sampling Line (Female Luer)	25
CO₂ Masks		
0562	Adult CO ₂ Mask w/ Swivel Elbow, .050 ID GSL	.50
0563	Adult Mask w/ CO ₂ Button	.50
Gas Sampling Accessories		
0535	GSL 10' PVC, .050" ID X /108" OD, Male luer, Swivel Luer Elbow w/Spin Lock Luer Cap	.50
8244	96" Adult expandable circuit, med. concentration Mask w/CO ₂ luer, B/V Filter	.20
Flex Extenders		
S6410P	- 4" Flex Tubing, 15mm OD x 15mm ID	50/cs
S6411P	- 6" Flex Tubing, 15mm OD x 15mm ID	50/cs
S6413P	- 6" Flex Tubing, 15mm OD x 22/15mm ID	50/cs
S6414P	- 4" 15mm OD x 22/15mm ID	50/cs
S6419P	- 6" Flex Tubing, 15mm x 22mm x 15mm	50/cs
S6420P	- 6" Flex Tubing, 15mm OD x 22/15 ID w/Capped Gas Sampling Port	50/cs
7-5132	- 3 1/2" Flex w/connector	50/cs

Westmed
 When Your Care is Critical

Westmed, Inc.

5580 S. Nogales Hwy.,
 Tucson, AZ 85706

Phone: 800-975-7987

Fax: 520-294-6061

www.westmedinc.com