

Hauge™ Airway Conscious Sedation Kit



Way Beyond an Airway...

Westmed

When Your Care is Critical

Hauge™ Airway for Total Airway Management

CONSCIOUS SEDATION KIT

The Hauge Airway is designed for simplicity and ease of use for most conscious sedation cases. The anatomically correct profile provides an immediate patent airway without stimulating a gag reflex. The Hauge Airway provides enough width to keep the tongue anterior in the oropharynx to ensure airway patency. Multiple applications for the disposable airway include:

- Awake Intubations
- Outpatient procedures
- Surgical MAC cases
- Upper Endoscopy procedures
- Bronchoscopy procedures
- Emergency room care
- For use with Esophageal Dilators
- Sleep Apnea Patients

Benefits for the provider:

- Immediate, established airway
- Universal O₂ and ETCO₂ ports
- Capability to deliver high flow O₂ orally to adult mouth breathers
- Large orifice accommodates a variety of scopes, suctioning
- Pre-procedure placement, with patient comfort in mind
- Device design accommodates elastic strap to ensure airway placement and stabilization
- Access to the oropharynx
- Bite block protection

Ordering Information:

- 6584 Large Hauge Airway 10/bx
- 6585 Medium Hauge Airway 10/bx
- 6586 Large Airway CSK 10/bx
- 6587 Medium Airway CSK 10/bx
- 6595 Hauge II Airway CSK 10/bx



Can accommodate:

- Suction catheter
- Fiber-optic scope with an endotracheal tube
- Upper Endoscopy scope
- TEE scope

Additional Uses:

- Nebulizer treatment – Corrugated tubing is connected to the large center opening
- Airway patency after general anesthesia

CSK includes:

- Hauge Airway (medium or large)
- CO₂/O₂ supply line
- Head strap to secure airway
- Tongue blade
- 2" Pigtail with Female luer

Westmed

When Your Care is Critical

Westmed, Inc.

5580 S. Nogales Hwy • Tucson, AZ 85706
Phone : 800/975-7987 • Fax : 520/294-6061
www.westmedinc.com • sales@westmedinc.com
Label PN : 74932-02