

# Uni-Lim™ & Pedi-Lim™

ADULT AND PEDIATRIC CO-AXIAL CIRCLE ANESTHESIA BREATHING CIRCUITS

The  
Optimal  
Choice



**Westmed**

*When Your Care is Critical*

# Uni-Lim™ & Pedi-Lim™

ADULT AND PEDIATRIC CO-AXIAL CIRCLE ANESTHESIA BREATHING CIRCUITS

The Westmed Uni-Lim™ is the optimal choice for your Co-Axial Circle System needs. Westmed offers the right product at an affordable price.

- Designed to function as a standard circle anesthesia breathing system.
- Features a streamlined tube-within-a-tube design which eliminates twisting between the inspiratory and expiratory limbs.
- Eliminates bulky connectors at the patient interface.
- Provides natural heat and moisture conservation properties to that of a standard circle circuit.
- Engineered with a secure-fit inspiratory limb to prevent patient end disconnects. No floating internal tubes.
- ISO-standard connections for secure, leak-resistant attachments.
- Available in convenient Latex-Free Configurations.
- Available in custom configurations to meet your clinical needs.
- Complete anesthesia kits available.

## Many Anesthesia Kit Options to Choose From


- Plain or ported wyes
- HME's, HMEF's, bacterial viral filters, and HEPA filters
- Standard and swivel gas sampling ports, lines and filters
- Internal and external gas sampling configurations
- Multiple anesthesia face mask options
- Multiple styles of anesthesia breathing bags
- Guedel and Berman airways
- ECG leads
- Flex extenders with and without gas sampling
- Many other options available

### Adult Breathing Circuits

Part Number	Description	Case Quantity
<b>Uni-Lim 40" Adult Co-Axial Circle Anesthesia Breathing Circuit</b>		
6651	GSE, 1 Filter, 3L Latex Bag . . . . .	20
6652	GSE, 1 Filter, 3L <b>Latex Free</b> Bag . . . . .	20
6655	GSE, 3L Latex Bag, Angled HMEF. . . . .	20
6656	GSE, 3L <b>Latex Free</b> Bag, Angled HMEF . . . . .	20
<b>Uni-Lim 60" Adult Co-Axial Circle Anesthesia Breathing Circuit</b>		
6653	GSE, 1 Filter, 3L Latex Bag . . . . .	20
6654	GSE, 1 Filter, 3L <b>Latex Free</b> Bag . . . . .	20

### Pediatric Breathing Circuits

<b>Pedi-Lim 40" Pediatric Co-Axial Circle Anesthesia Breathing Circuit</b>		
8121	GSE, 1 Filter, 1L <b>Latex Free</b> Bag . . . . .	20
<b>Pedi-Lim 60" Pediatric Co-Axial Circle Anesthesia Breathing Circuit</b>		
8122	GSE, 1 Filter, 1L <b>Latex Free</b> Bag . . . . .	20



Ask your local representative about other standard and custom configurations that are available to meet your breathing circuit needs.

Please log onto [www.westmedinc.com](http://www.westmedinc.com) for a complete listing of Westmed anesthesia and respiratory products.

© 2007 Westmed, Inc. All rights reserved. PN 74792 (Rev.01)

**Westmed**  
When Your Care is Critical

**Westmed, Inc.**  
5580 S. Nogales Hwy.,  
Tucson, AZ 85706  
Phone: 800-975-7987  
Fax: 520-294-6061  
[www.westmedinc.com](http://www.westmedinc.com)